

Product Datasheet

SFRS12 antibody (orb317845)

Description	SFRS12 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human SFRS12/p18SRP (21-120/155 aa)
Target	SFRS12/p18SRP
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18 kDa
Uniprot ID	Q8N9Q2
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.