


Product Datasheet

SLC39A13 antibody (orb317967)


Description	SLC39A13 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human SLC39A13 (1-100/371 aa)
Target	SLC39A13
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q96H72
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Sample: Small intestine
(Mouse) Lysate
a...



Sample:
U87MG(Human)
Cell Lysate at
30 u...



Sample:
K562(Human)
Cell Lysate at
30 ug...