



Product Datasheet

TSPAN6 antibody (orb317988)


Description	TSPAN6 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human TSPAN6 (1-100/245 aa)
Target	TSPAN6
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28 kDa
Uniprot ID	O43657
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse hepatoma cells



Immunohistochemical staining of mouse hepatoma cells



Flow cytometric analysis of HepG2 cells