

Product Datasheet

VGLUT3 antibody (orb317998)

Description	VGLUT3 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human VGLUT3/SLC17A8 (1-100/589 aa)
Target	VGLUT3/SLC17A8
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65 kDa
Uniprot ID	Q8NDX2
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.