

## Product Datasheet

### ZIP Kinase (phospho-Thr265) antibody (orb318041)

|                            |                                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | ZIP Kinase (phospho-Thr265) antibody                                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                                                        |
| <b>Reactivity</b>          | Mouse, Rat                                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                                  |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                     |
| <b>Immunogen</b>           | KLH conjugated synthesised phosphopeptide derived from human ZIP Kinase around the phosphorylation site of Thr265 (RM(p-T)IA) |
| <b>Target</b>              | phospho-ZIP Kinase (Thr265)                                                                                                   |
| <b>Preservatives</b>       | Preservative: 15mM Sodium Azide,<br>Constituents: 1% rAlbumin, 0.01M PBS,<br>pH 7.4                                           |
| <b>Form/Appearance</b>     | Liquid                                                                                                                        |
| <b>Concentration</b>       | 1mg/ml                                                                                                                        |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                                              |
| <b>Note</b>                | For research use only                                                                                                         |
| <b>Isotype</b>             | IgG                                                                                                                           |
| <b>Clonality</b>           | Polyclonal                                                                                                                    |
| <b>MW</b>                  | 53 kDa                                                                                                                        |
| <b>Uniprot ID</b>          | <b>043293</b>                                                                                                                 |
| <b>Dilution Range</b>      | WB=1:500-2000, IHC-P=1:100-500                                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                               |