

Product Datasheet

TBC1D3C antibody (FITC) (orb31806)

Description	TBC1D3C antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human TBC1 domain family member 3C protein (1-260AA)
Target	TBC1D3C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6IPX1
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.