

Product Datasheet

Leptin Receptor antibody (orb318718)

Description

Leptin Receptor antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Leptin Receptor. Please contact us for the exact immunogen sequence. The peptide is available as orb374915.

Target

Leptin Receptor

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

132 kDa

Uniprot ID

[P48356](#), [Q62959](#), [P48357](#)

NCBI

[001003679](#)

Entrez



[23484](#)

Dilution Range


IHC-P: 1:100-800, WB: 1:100-800

Expiration Date


12 months from date of receipt.



IHC-P image of pig small intestines tiss...



Immunohistochemical staining of mouse br...



Western blot analysis of mouse liver (la...