


Product Datasheet


Sialoadhesin antibody (orb318719)

Description	Sialoadhesin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Sialoadhesin. Please contact us for the exact immunogen sequence. The peptide is available as orb374916.
Target	Sialoadhesin
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	Q62230 , Q9BZZ2
NCBI	075556 , 023068
Entrez	6614
Dilution Range	IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.




IHC-P image of rat bladder tissue using

...



Immunohistochemical staining of paraffin...



IHC-P staining of mouse lymph node tissue...