

Product Datasheet

iNOS antibody (orb318724)

Description

iNOS antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human iNOS. Please contact us for the exact immunogen sequence. The peptide is available as orb12914.

Target

iNOS

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

131 kDa

Uniprot ID

[P35228](#), [Q06518](#), [P29477](#)

NCBI

[31](#), [41](#), [000625](#), [000616](#)

Entrez


[4843](#)

Dilution Range


IHC-P: 1:100-600, WB: 1:100-800

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of rat lung...



Western blot analysis of mouse lung (lan...



IHC-P image of mouse brain tissue using ...