


Product Datasheet

CHRNA7 antibody (orb318729)

Description	Goat polyclonal antibody to CHRNA7
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	CHRNA7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Rat Cerebral Cortex and Rat Hippocampus shows membranous staining in neural cells. Recommended concentration: 5-10µg/ml. Paraffin embedded Human Brain (Cortex). Recommended concentration: 3-6µg/ml.
Protein Sequence	KRPGEDKVRPACQHKQ
Clonality	Polyclonal
MW	56.5
Entrez	1139
Dilution Range	ELISA: 1:128000, IHC-P: 3-6 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Rat Cereb...