



Product Datasheet

BDNF antibody (orb318730)


Description	Goat polyclonal antibody to BDNF
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	BDNF
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Immunofluorescence: Expression of the protein seen in the nucleus and cytoplasm of Glioblastoma U251 cells. Recommended concentration: 10µg/ml.
Protein Sequence	ETKCNPMGYTKE
Clonality	Polyclonal
MW	27.8; 29.8
Entrez	627
Dilution Range	ELISA: 1:2000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Rat Cerebellum



Immunofluorescence analysis of Mouse sciatic nerve



Immunohistochemical staining of Rat Cerebellum