


Product Datasheet

CYP26A1 antibody (orb318731)

Description	Goat polyclonal antibody to CYP26A1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC
Target	CYP26A1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Brain (Cortex and Cerebellum). Recommended concentration: 3.75µg/ml.. Notes: Immunocytochemistry: ER-staining of PFA-fixed and TX100-treated HeLa cells. Recommended concentration: 2-5ug/ml.
Protein Sequence	NLPARFTHFHGE
Clonality	Polyclonal
MW	56.2; 48.6
Entrez	1592
Dilution Range	ELISA: 1:16000, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...