


Product Datasheet

CYP2B6 antibody (orb318732)

Description	Goat polyclonal antibody to CYP2B6
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CYP2B6 (aa247-259)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 75kDa band observed in lysates of cell line HepG2 and K562 (calculated MW of 56.3kDa according to NP_000758.1). Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	HSVEKHRETLDPS
Clonality	Polyclonal
MW	56.3
Entrez	1555
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 (Lane 1) ...