

Product Datasheet

CSNK2B antibody (orb318735)

Description

Goat polyclonal antibody to CSNK2B

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Target

CSNK2B

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

This antibody is expected to recognize both reported isoforms (NP_001311.3; NP_001269314.1). ELISA: antibody detection limit dilution 1:128000. WB: Approx 26kDa band observed in lysates of cell lines NIH3T3, HeLa, K562, MCF7 and HepG2 (calculated MW of 24.9kDa according to NP_001311.3). Approx 25kDa band observed in Mouse Testis, 26kDa in Rat Testis and 28kDa in Pig Testis lysates (calculated MW of 24.9kDa according to Mouse NP_001290374.1, Rat NP_001030315.1 and pig NP_001138849.1). The observed molecular weights correspond to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.3-1 µg/ml.

Protein Sequence

TGLNEQVPYHRQ

Clonality

Polyclonal

MW

24.9

Entrez


1460

Dilution Range


ELISA: 1:128000, WB: 0.3-1 µg/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of mouse (Lane1), ...



Western blot analysis of Mouse, Rat and ...