



Product Datasheet

GNA12 antibody (orb318738)

Description	Goat polyclonal antibody to GNA12
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	GNA12
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize all reported isoforms (NP_031379.2; NP_001269369.1; NP_001269370.1; NP_001280021.1). ELISA: antibody detection limit dilution 1:128000. WB: Approx 37kDa band observed in Human Placenta lysates and approx 42kDa in Mouse Heart lysates (calculated MW of 37.6kDa according to Human NP_001269370.1 and 44.1kDa according to Mouse NP_034432.1). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	KKHFPDFRGDPHR
Clonality	Polyclonal
MW	44.3; 35.2; 37.6; 42.3
Entrez	2768
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Placenta ...



Western blot analysis of Mouse Heart lys...