

Product Datasheet

BOD1L antibody (orb318739)

Description	Goat polyclonal antibody to BOD1L
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	BOD1L
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 330kDa according to NP_683692.2. Therefore we cannot recommend an optimal concentration and the antibody is an aspiring p
Protein Sequence	RTRGQQRVEEAP
Clonality	Polyclonal
MW	330
Entrez	259282
Dilution Range	ELISA: 1:128000
Expiration Date	12 months from date of receipt.