


Product Datasheet

CD87 antibody (orb318741)

Description	Goat polyclonal antibody to CD87
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	U-PAR / CD87
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 60kDa band observed in lysates of cell line U251 (calculated MW of 37.0kDa according to NP_002650.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 1-3µg/ml.
Protein Sequence	RTGLKITSLTE
Clonality	Polyclonal
MW	37.0; 31.3; 32.0
Entrez	5329
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of U251 lysate using...