


Product Datasheet

HSPE1 antibody (orb318743)

Description	Goat polyclonal antibody to HSPE1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CPN10 / HSPE1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 10kDa band observed in lysates of cell line HepG2 (calculated MW of 10.9kDa according to NP_002148.1). Recommended concentration: 0.01-0.03µg/ml.. IHC: Paraffin embedded Human Liver. Recommended concentration: 5µg/ml.
Protein Sequence	PEKSQGKVLQAT
Clonality	Polyclonal
MW	10.9
Entrez	3336 , 100529241
Dilution Range	ELISA: 1:128000, WB: 0.01-0.03 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate us...