



Product Datasheet

HLA-B antibody (orb318746)

Description	Goat polyclonal antibody to HLA-B
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	HLA-B (aa132-143)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 40+45kDa bands observed in lysates of cell line Jurkat (calculated MW of 40.5kDa according to NP_005505.2). The observed 45kDa band corresponds to the glycosylated form. Recommended concentration: 0.1-0.3µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Tonsil. Recommended concentration: 2.5µg/ml.
Protein Sequence	RLLRGHDQYAYD
Clonality	Polyclonal
MW	40.5
Entrez	3106
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...



Immunohistochemical staining of Human To...