


Product Datasheet

IL13RA2 antibody (orb318751)

Description	Goat polyclonal antibody to IL13RA2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	IL13RA2 (aa96-107)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 50+60kDa bands observed in lysates of cell line HeLa (calculated MW of 44.2kDa according to NP_000631.1). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Kioi et al, FASEB J. 2006 Nov;20(13):23
Protein Sequence	DLNKGIEAKIHT
Clonality	Polyclonal
MW	44.2
Entrez	3598
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using IL13RA2 antibody (orb318751).