



## Product Datasheet

### ZRANB2 antibody (orb318769)


|                            |                                                                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ZNF265                                                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                                                                                |
| <b>Reactivity</b>          | Human, Mouse, Porcine, Primate, Rabbit, Rat                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                                          |
| <b>Tested Applications</b> | IF, IH, WB                                                                                                                            |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ZNF265. The exact sequence is proprietary. |
| <b>Target</b>              | ZRANB2                                                                                                                                |
| <b>Preservatives</b>       | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                             |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                             |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                      |
| <b>Note</b>                | For research use only                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                            |
| <b>Source</b>              | Rabbit                                                                                                                                |
| <b>Uniprot ID</b>          | <a href="#">O35986</a> , <a href="#">Q9R020</a> , <a href="#">O95218</a>                                                              |
| <b>Entrez</b>              | <a href="#">53861</a> , <a href="#">58821</a> , <a href="#">9406</a>                                                                  |
| <b>Dilution Range</b>      | WB: 1:500-1:1000, IHC-P: 1:100-1:200, IF/ICC: 1:100-1:500                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                       |



Western blot analysis of ZNF265 expressi...



Immunohistochemical analysis of ZNF265 S...



Immunofluorescent analysis of ZNF265 sta...