



Product Datasheet

H2AFX (phospho-S139) antibody (orb318857)


Description	Rabbit polyclonal antibody to Histone H2A.X (Phospho-S139)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S139 of human Histone H2A.X protein. The exact sequence is proprietary.
Target	H2AFX
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P16104, P27661
Entrez	15270, 3014
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200, IP: 1:10-1:100, FACS: 1:100-1:200M
Expiration Date	12 months from date of receipt.



Western blot analysis of Histone H2A.X (...)



Immunohistochemical analysis of Histone ...



Immunofluorescent analysis of Histone H2...