



Product Datasheet

TUBA1A antibody (orb318873)


Description	Rabbit polyclonal antibody to Alpha-tubulin
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human Alpha-tubulin. The exact sequence is proprietary.
Target	TUBA1A
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P68369 , Q71U36 , P68370 , Q6P9V9 , P68363 , P05213
Entrez	500929 , 22142 , 64158 , 22143 , 10376 , 7846
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:100-1:200, IF/ICC: 1:100-1:500, IP: 1:10-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of Alpha-tubulin e...



Immunohistochemical analysis of Alpha-tu...



Immunofluorescent analysis of Alpha-tubu...