

Product Datasheet

EGFR (phospho-S1026) peptide (orb319106)

Description The peptide is used to anti EGFR (phospho-S1026) antibody reactivity.

Reactivity Human

Tested Applications Blocking

Target EGFR

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes This peptide was used as immunogen to produce orb318765 (EGFR (phospho-S1026) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID **P00533**

Entrez **1956**

Expiration Date 12 months from date of receipt.