

Product Datasheet

ARHGAP24 antibody (orb31911)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ARHGAP24
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	84258
Uniprot ID	Q8N264
NCBI	NP_001036134.1 , NP_001020787.2 , NP_112595.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

