

Product Datasheet

FOS (phospho-S32) peptide (orb319203)

Description	The peptide is used to anti FOS (phospho-S32) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	FOS
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	<p>This peptide was used as immunogen to produce orb318862 (FOS (phospho-S32) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.</p>
Purity	> 85%
Source	Synthetic
Uniprot ID	P01100, P12841, P01101
Entrez	14281, 314322, 2353
Expiration Date	12 months from date of receipt.