

Product Datasheet

TUBA1A peptide (orb319214)

Description	The peptide is used to anti TUBA1A antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	TUBA1A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb318873 (TUBA1A antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05213 , P68370 , P68363 , Q71U36 , P68369 , Q6P9V9
Entrez	64158 , 10376 , 500929 , 22142 , 22143 , 7846
Expiration Date	12 months from date of receipt.