

Product Datasheet

RHOBTB3 peptide (orb319227)

Description The peptide is used to anti RHOBTB3 antibody reactivity.

Reactivity Human, Mouse

Tested Applications Blocking

Target RHOBTB3

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes This peptide was used as immunogen to produce orb318886 (RHOBTB3 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID [Q9CTN4](#), [O94955](#)

Entrez [22836](#), [73296](#)

Expiration Date 12 months from date of receipt.