

Product Datasheet

EP300 AcK1558/K1560) peptide (orb319232)

Description	The peptide is used to anti EP300 AcK1558/K1560) antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	EP300
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb318891 (EP300 AcK1558/K1560) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	B2RWS6, Q09472
Entrez	2033, 328572
Expiration Date	12 months from date of receipt.