

Product Datasheet

LMO3 peptide (orb319289)

Description	The peptide is used to anti LMO3 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	LMO3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb318948 (LMO3 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8TAP4, Q8BZL8, Q99MB5
Entrez	55885, 109593
Expiration Date	12 months from date of receipt.