

Product Datasheet

ASF1B peptide (orb319318)

Description	The peptide is used to anti ASF1B antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	ASF1B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb318977 (ASF1B antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9NVP2
Entrez	55723
Expiration Date	12 months from date of receipt.