

Product Datasheet

SVOP peptide (orb319385)

| | |
|----------------------------|---|
| Description | The peptide is used to anti SVOP antibody reactivity. |
| Reactivity | Human, Mouse, Rat |
| Tested Applications | Blocking |
| Target | SVOP |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | This peptide was used as immunogen to produce orb319044 (SVOP antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | Q8BFT9 , Q8N4V2 , Q9Z2I7 |
| Entrez | 55530 , 171442 , 68666 |
| Expiration Date | 12 months from date of receipt. |