

Product Datasheet

UBE3B peptide (orb319392)

Description	The peptide is used to anti UBE3B antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	UBE3B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	<p>This peptide was used as immunogen to produce orb319051 (UBE3B antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.</p>
Purity	> 85%
Source	Synthetic
Uniprot ID	Q7Z3V4
Entrez	89910
Expiration Date	12 months from date of receipt.