

## Product Datasheet

### SERPINA3 peptide (orb319411)

**Description** The peptide is used to anti SERPINA3 antibody reactivity.

**Reactivity** Human

**Tested Applications** Blocking

**Target** SERPINA3

**Form/Appearance** Lyophilized powder

**Storage** Shipped at 4°C. Store at -20°C for one year.

**Note** For research use only

**Application notes** This peptide was used as immunogen to produce orb319070 (SERPINA3 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

**Purity** > 85%

**Source** Synthetic

**Uniprot ID** **P01011**

**Entrez** **12**

**Expiration Date** 12 months from date of receipt.