

## Product Datasheet

### TMEM145 peptide (orb319430)

**Description** The peptide is used to anti TMEM145 antibody reactivity.

**Reactivity** Human, Mouse

**Tested Applications** Blocking

**Target** TMEM145

**Form/Appearance** Lyophilized powder

**Storage** Shipped at 4°C. Store at -20°C for one year.

**Note** For research use only

**Application notes** This peptide was used as immunogen to produce orb319089 (TMEM145 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

**Purity** > 85%

**Source** Synthetic

**Uniprot ID** [Q8NBT3](#), [Q8C4U2](#)

**Entrez** [330485](#), [284339](#)

**Expiration Date** 12 months from date of receipt.