

Product Datasheet

IL10 antibody (orb319439)

Description

IL10 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications
Immunogen

KLH conjugated synthetic peptide derived from human IL10. Please contact us for the exact immunogen sequence. The peptide is available as orb374925.

Target

IL10. Recognizes endogenous levels of IL10 protein

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

21 kDa

Uniprot ID

[P22301](#), [P29456](#), [P18893](#)

NCBI

[000572](#), [000563](#)

Entrez



[3586](#)

Dilution Range


IHC-P: 1:100-500 (based on 0.5 mg/ml), WB: 1:100-500

Expiration Date


12 months from date of receipt.



IHC-P staining of pig small intestines t...



Immunohistochemical staining of pig smal...



IHC-P image of rat lymph node tissue usi...