

Product Datasheet

CD105 antibody (orb319440)

Description

CD105 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human CD105. Please contact us for the exact immunogen sequence. The peptide is available as orb374926.

Target

CD105. Recognizes endogenous levels of CD105 protein

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

71 kDa

Uniprot ID

[P17813](#), [Q63961](#)

NCBI

[000118](#), [31](#), [000109](#), [11](#)

Entrez


[2022](#)

Dilution Range

IHC-P: 1:100-800, WB: 1:100-500


Expiration Date

12 months from date of receipt.




WB analysis of rat lymph node (lane 1),

...



IHC-P image of rat lymph node tissue using



Immunohistochemical staining of mouse brain tissue