

Product Datasheet

CD24 antibody (orb319443)

Description

CD24 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical staining of paraffin...

Immunogen

KLH conjugated synthetic peptide derived from human CD24. Please contact us for the exact immunogen sequence. The peptide is available as orb374929.

Target

CD24. Recognizes endogenous levels of CD24 protein

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

IHC-P image of rat lymph node tissue usi...

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

Immunohistochemical staining of paraffin...

MW

11 kDa

Uniprot ID
P25063
Entrez
100133941
Dilution Range

IHC-P: 1:100-500 (based on 0.5 mg/ml)

Expiration Date

12 months from date of receipt.