


Product Datasheet

CD150 antibody (orb319446)

Description	CD150 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD150. Please contact us for the exact immunogen sequence. The peptide is available as orb374932.
Target	CD150. Recognizes endogenous levels of CD150 protein
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	37 kDa
Uniprot ID	Q13291, Q9QUM4
Entrez	6504
Dilution Range	WB: 1:100-500 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Western blot analysis of rat colon (lane...)