



Product Datasheet

MCU antibody (orb319447)


Description	MCU antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human MCU. Please contact us for the exact immunogen sequence. The peptide is available as orb374933.
Target	MCU. Recognizes endogenous levels of MCU protein
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Dilution Range	IHC-P: 1:100-800 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



IHC-P staining of rat stomach tissue using MCU antibody.



Immunohistochemical staining of rat testis using MCU antibody.



Immunohistochemical staining of paraffin-embedded tissue using MCU antibody.