

Product Datasheet

Alpha Synuclein (ser129) Antibody (orb319462)

Description

Our Anti-Alpha Synuclein (Ser129) rabbit

antibody from Biorbyt is...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser129 of rat alpha-synuclein, conjugated to keyhole limpet hemocyanin (KLH).

Target

Alpha Synuclein (ser129)

Preservatives

10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml rAlbumin and 50% glycerol.

Form/Appearance

Antigen Affinity Purified from Pooled Serum

Storage

Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

15

Uniprot ID

P37377

Dilution Range

WB: 1:1000, IHC: 1:200

Expiration Date

12 months from date of receipt.