

Product Datasheet

COBRA1 (NELF-B) Antibody (orb319483)

Description	Our Anti-COBRA1 (NELF-B) primary antibody from Biorbyt is mouse monoclonal. It detects human, mouse,...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant full length human COBRA1 purified from E. coli.
Target	COBRA1 (NELF-B)
Preservatives	PBS + 0.05% NaN ₃
Form/Appearance	Protein G Purified
Storage	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	62
Uniprot ID	Q8WX92
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.