

Product Datasheet

COXIV (Ser58) Antibody (orb319489)

Description	Our Anti-COXIV (Ser58) rabbit polyclonal phosphospecific primary antibody from Biorbyt is produced i...
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser58 of mouse COXIV, conjugated to keyhole limpet hemocyanin (KLH).
Target	COXIV (Ser58)
Preservatives	10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml rAlbumin and 50% glycerol.
Form/Appearance	Antigen Affinity Purified from Pooled Serum
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17
Uniprot ID	P19783
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.