

Product Datasheet

DARPP-32 (Thr75) Antibody (orb319493)

Description	Our Anti-DARPP-32 (Thr75) phosphospecific primary antibody from Biorbyt is rabbit polyclonal. It det...
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phospho-peptide corresponding to amino acid residues surrounding Thr75 of rat DARPP-32 conjugated to keyhole limpet hemocyanin (KLH).
Target	DARPP-32 (Thr75)
Preservatives	10 mM HEPES (pH 7.5), 150 mM NaCl, 100 µg per ml rAlbumin and 50% glycerol.
Form/Appearance	Antigen Affinity Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32
Uniprot ID	Q6J4I0
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.