



Product Datasheet

Doublecortin Antibody (orb319504)


Description	Our Anti-Doublecortin primary antibody from Biorbyt is mouse monoclonal. It detects bovine, human, m...
Species/Host	Mouse
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Full length recombinant human protein expressed in and purified from E. coli.
Target	Doublecortin
Preservatives	PBS + 50% glycerol and 5 mm NaN3
Form/Appearance	Protein G Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG2A
Clonality	Monoclonal
MW	35, 45
Uniprot ID	O43602
Dilution Range	WB: 1:5000, IHC: 1:1000, ICC: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of rat neuro...



Immunofluorescence analysis of rat corti...



Western blot analysis of rat brain lysat...