

Product Datasheet

GABAB Receptor (Ser923), R1-Subunit Antibody (orb319543)

Description

Our Anti-GABAB Receptor (Ser923),

phosphospecific primary antibody
from...

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Synthetic phospho-peptide
corresponding to amino acid residues
surrounding Ser923 of the R1-subunit
of rat GABAB receptor, conjugated to
keyhole limpet hemocyanin (KLH).

Target

GABAB Receptor (Ser923), R1-
Subunit

Preservatives

10 mM HEPES (pH 7.5), 150 mM
NaCl, 100 µg per ml rAlbumin and
50% glycerol.

Form/Appearance

Antigen Affinity Purified from Pooled
Serum

Storage

Storage at -20°C is recommended, as
aliquots may be taken without
freeze/thawing due to presence of
50% glycerol. Stable for at least 1
year at -20°C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

102

Uniprot ID

Q9Z0U4

Dilution Range

WB: 1:1000

Expiration Date

12 months from date of receipt.