

Product Datasheet

Glycogen synthase kinase 3 B (GSK3B) (Ser9) Antibody (orb319559)

Description

Our Anti-Glycogen synthase kinase

	polyclonal phosphospecific primary antib...
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser9 of rat GSK3 β , conjugated to keyhole limpet hemocyanin (KLH).
Target	Glycogen synthase kinase 3 B (GSK3B) (Ser9)
Preservatives	10 mM HEPES (pH 7.5), 150 mM NaCl, 100 μ g per ml rAlbumin and 50% glycerol.
Form/Appearance	Antigen Affinity Purified from Pooled Serum
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46
Uniprot ID	P18266
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.