

Product Datasheet

KCNC2 antibody (orb31956)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KCNC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	70226
Uniprot ID	Q96PR1
NCBI	NP_631875.1 , NP_001247428.1 , NP_001247427.1 , NP_715624.1 , NP_631874.1 , NP_001247426.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

