



Product Datasheet

LAMP1 Antibody (orb319565)


Description	Our Anti-LAMP1 primary antibody from Biorbyt is mouse monoclonal. It detects human LAMP1 and is Prot...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Recombinant full length human LAMP1 expressed in and purified from E. coli.
Target	LAMP1
Preservatives	PBS + 10 mM NaN ₃ .
Form/Appearance	Protein G Purified
Storage	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	100
Uniprot ID	P11279
Dilution Range	WB: 1:5000, IHC: 1:2000, ICC: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunostaining analysis of HeLa cells lysates.



Immunolabeling analysis of HeLa cells lysates.



Western blot analysis of HeLa cell lysates.