

Product Datasheet

Microtubule Associated Protein 2 (MAP2) Antibody (orb319566)

Description

Our Anti-Microtubule Associated

from Biorbyt is chicken polyclonal...

Species/Host

Gallus

Reactivity

Bovine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC, WB

Applications

Immunogen

Recombinant bovine MAP2 protein expressed in and purified from E. Coli.

Target

Microtubule Associated Protein 2 (MAP2)

Preservatives

PBS + 10 mM NaN₃.

Form/Appearance

Total IgY fraction

Storage

Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20° C or -80° C, while keeping a working aliquot stored at 4° C for short term. Avoid freeze/thaw cycles. Stable for at least 1 year.

Note

For research use only

Isotype

IgY

Clonality

Polyclonal

MW

280

Uniprot ID

P11137

Dilution Range

WB: 1:20000, IHC: 1:2500-1:10,000, ICC: 1:1000-1:20,000

Expiration Date

12 months from date of receipt.